

Allen.
Hug-U-Vac.
Steep Trend Positioners



Ideal for Use with Robotic Procedures

The Allen Hug-U-Vac Steep Trend Positioners safely and comfortably hold patients in the steep Trendelenburg. The Hug-U-Vac Positioners have integral quick-release straps that enable it to attach tightly to the rails of standard operating tables.

In addition, the Allen Hug-U-Vac Positioner's unique pliable design enables it to inflate and deflate evenly, to fit the patient's body. The Large Hug-U-Vac Steep Trend Positioner has been designed with three chambers to accommodate larger patients.

Ideal for use with:

- Gynecology Procedures
- Colorectal Procedures
- Urology Procedures
- Cholecystectomies
- Robotic Surgical Procedures

At a glance

- Designed for easy forearm access to IV lines
- Patient weight capacity: 500 lbs (227 kg)
- Enables secure positioning in robotic procedures
- Integral straps secure to OR table rails
- Saves positioning time
- Valve location designed for quick access to anesthesia
- Lightweight for easy transport, set-up, and storage
- No continuous suction required
- Durable vacuum valve with on/off switch and lock
- Safe for use with warming blankets
- Disposable slipcovers are waterproof for easy clean-up

Allen.
Hug-U-Vac.
Steep Trend Positioners

About the Product

The Allen Hug-U-Vac Steep Trend Positioners are used for positioning patients in the steep Trendelenburg position for surgical procedures. They have a patient weight capacity of 500 lbs (227 kg). The Hug-U-Vac Positioners fit standard OR tables, by securing its three straps to the rails of the OR table.

Allen's Hug-U-Vac Steep Trend Positioners are filled with soft microbeads that enable it to inflate and deflate evenly, when suction is applied. No continuous suction is required. After positioning the patient, deflate the device and its pliable shape conforms to the contours of the patient's body. Allen Hug-U-Vac Steep Trend disposables cover the device completely, and enable easy clean-up to expedite case turnover. The Steep Trend Positioners should be used with Allen Steep Trend Shoulder Supports to ensure additional safety in the steep Trendelenburg position.



Allen Hug-U-Vac Steep Trend Positioner



Allen Hug-U-Vac Large Steep Trend Positioner



Steep Trend Positioner Disposable Covers



Integral Straps



Safety Vacuum Valve



Allen Steep Trend Shoulder Supports

Technical Specifications

Patient Weight Capacity 500 lbs (227 kg)

Allen Hug-U-Vac Steep Trend Positioner

Device Dimensions..... 46 in x 42 in (116.8 cm x 106.7 cm)

Device Weight..... 8 lbs (3.6 kg)

Allen Hug-U-Vac Large Steep Trend Positioner

Device Dimensions..... 51 in x 54 in (129.5 cm x 137.1 cm)

Device Weight..... 10 lbs (4.5 kg)

#A-60010 Steep Trend Positioners Package	\$3,749
Includes Hug-U-Vac Steep Trend Positioner and Hug-U-Vac® Large Steep Trend Positioner	
#A-60001 Allen Hug-U-Vac Steep Trend Positioner	\$1,749
#A-60005 Allen Hug-U-Vac Large Steep Trend Positioner	\$2,499
#A-60101 Allen Hug-U-Vac Steep Trend Disposable Covers (Case of 12) (shown)	\$220
#A-60105 Allen Hug-U-Vac Large Steep Trend Disposable Covers (Case of 12)	\$270
#A-63600 Allen Steep Trend Shoulder Supports (works with Steep Trend and Large Steep Trend Positioners)	\$999

Allen® is a registered trademark of Allen Medical Systems. Hug-U-Vac® is a registered trademark of Hug-U-Vac Surgical Positioning Systems, Inc.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2012 Allen Medical, A Hill-Rom Company.
 ALL RIGHTS RESERVED.



A **Hill-Rom** Company

D-770599-A3 PENDING

USA 800-433-5774
 www.allenmedical.com