


# S U R G I C A L   T A B L E S

## **E·** mayfield® Infant Cast Table & Traction Conversion Kit



- Designed to be versatile; with or without traction
- Full body casting, hip and shoulder spicas
- Compact and light weight

# mayfield® Infant Cast Table & Traction Conversion Kit

## It's adaptable

With traction or without. For Infants or children. The Infant Cast Table becomes just what you want it to be, with all the accessories you need for a multitude of applications.

## The Infant Cast Table With Traction

It's designed and sized specifically for the application of casts on infants and children. The traction model is ideal for hip spicas. It uses a simple expansion principal and saddle changes that allow easy conversion from infant to child set-up.

## Hip Spicas, Shoulder Spicas

It's designed to be versatile. Use the Infant Cast Table for full body casting and hip and shoulder spicas.

## For Infant or Child

An adjustable forward section accommodates infants and children. Child and infant saddles and perineal posts are interchangeable.

## Special Feature Head Section

Complete 180° rotation of the head section provide shoulder support or body support. An optional overhead frame is also available for the Infant Cast Table with Traction.

## More Versatility

The Infant Cast Table can be used on a stretcher, on an examining table, or in a crib. This choice of settings makes it easier on the patient- and easier on the staff.

You get a minimum of removable parts for a better contained, easier-to-handle unit. And when the procedure is over, the staff will find that the Infant Cast Table's light weight and compact design make it easy to store.


## Convert to Traction Later

Should you decide to purchase your table without traction, you can still convert to traction later. Just request the Traction Conversion Kit. Two leg spars and two traction units make it easy.

## Infant Cast Table

- 10005 Basic unit includes base table with head section and:
- 10209 Radiolucent, perineal post, 1.25" (3.2 cm)
- 10354 Cast saddle, child
- 22238 Infant body slab
- 22236 Infant saddle
- 22237 Metal perineal post, .875" (2.2 cm)
- 22239 Adapter sleeve for conversion kit


Shipping weight: 35 pounds (15.9 kg)



## Infant Cast Table with Traction


- 10006 Basic unit includes base table with head section, traction unit
- 10209 Radiolucent perineal post, 1.25" (3.2 cm)
- 10354 Cast saddle, child
- 10385 Traction conversion kit
- 22238 Infant body slab
- 22236 Infant saddle
- 22237 Metal perineal post, .875" (2.2 cm)
- 22239 Adapter sleeve for conversion kit

Shipping weight: 45 pounds (20.4 kg)



## Traction Conversion Kit

- 10385 Basic unit includes two each leg spars with traction unit
- Shipping weight: 25 pounds (11.4 kg)



 **schaerermayfield™**  
WHEREVER YOU OPERATE