

Allen.
Easy Catcher™
Disposable System



Easy to Use Fluid Collection System with an Integrated, Disposable Band

The Allen Easy Catcher Disposable System is easy to use and fully disposable. It mounts directly onto operating table rails and is disposed of after a surgical case. The Easy Catcher system contains a flexible band integrated into the drain bag, eliminating the need for separate bands or additional securing methods.

Ideal for use with cystoscopy and hysteroscopy procedures.

At a glance

- Integrated band and drain bag – eliminates the need for separate bands
- Pre-assembled system – eliminates the need to remove bag from the band
- Flexible band – enables unobstructed surgical site access
- Buttocks flap – ensures drainage directly into the catcher
- Integrated mesh trap – ensures proper fluid flow
- Removable hose adaptor – enables compatibility with multiple containment systems



A Hill-Rom Company

About the Product

The Allen Easy Catcher system is a fluid collection system for urology and gynecology procedures that attaches to the rails of an operating table. The drain bag includes a flap designed to slide under the patient to direct fluid into the bag, a mesh trap for proper drainage, and an extendable hose to attach to a fluid containment system.

Unlike other systems, the Easy Catcher system includes a disposable band integrated into the drain bag. The robust band design keeps the drain bag open for optimal fluid drainage and flexes in to allow unobstructed surgeon access.

Technical Specifications

Individually Pouched Disposable
 Fits Standard Operating Table Widths
 Latex-free

Drain hose:
 Expands to 5 ft (1.5 m)
 7/8" (22.2 mm) Diameter

Adaptor:
 1/4" (6.4 mm) Diameter



Allen Easy Catcher



Integral Band and Drain Bag



Mounts to Operating Table Rails



Unobstructed Surgeon Access



Drain Hose and Adaptor



Fully Disposable

#O-EC51 Easy Catcher Disposable System: STERILE (Case of 10) (US Rail) \$195

#O-EC52 Easy Catcher Disposable System: NON-STERILE (Case of 10) (US Rail) \$175

Allen® is a registered trademark of Allen Medical Systems.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2012 Allen Medical, A Hill-Rom Company.
 ALL RIGHTS RESERVED.

D-770609-A2



USA 800-433-5774
www.allenmedical.com