

Allen Shoulder Access System

Unlimited Surgical Site Access In The Sterile Field

The Allen Shoulder Access System is the smallest, lightweight arm positioner on the market – for any shoulder procedure in the Beach Chair position. It attaches over the drape and provides effortless positioning of the arm intra-operatively - at a fraction of the price of the Spider and the McConnell. Additionally, the disposable is more intuitive and competitively priced.



Ideal For Use With...

All shoulder and arthroscopic procedures including:

- Rotator Cuff Repair
- Bankart Repair
- SLAP Lesion Repair
- Capsular Release

About The System...

Spend your time on the procedure, not the positioning!

The Allen Shoulder Access System delivers nearly effortless arm positioning for shoulder surgery in the Beach Chair position. By indexing at the elbow, the device is able to gain complete positioning control of the shoulder joint. The revolutionary patent-pending design allows the entire device to be steam sterilized, and provides an unprecedented level of access to the surgical site. The device also has the capability of achieving a semi locked state, allowing the surgeon to simply move the arm without having to adjust the tension.

Key Features & Benefits

- Ease of Use: Effortless Positioning
“Just Move It”: Semi-locked state allows repositioning without tension adjustment
- Sterility: Connects over the drape, in the sterile field
- Compact storage and easy portability
- Single knob releases all joints
- Index at the elbow for complete control
- Clamp included
- The only fully autoclavable shoulder positioner on the market.
- Quick set-up
- Close access to the surgical site
- No Foot Pedal / No Bruised Shins
- No Crowding

Specifications

- Fits the Following Rails: US, EU, UK, DEN (specify when ordering)
- Device weighs only 6 lbs. (2.7 kg)
- Can support 500 lb. (227 kg) maximum patient weight

#A-92000	Shoulder Access System (US Rail)	\$4,995
#A-92000E	Shoulder Access System (EU Rail)	\$4,995
#A-92000UK	Shoulder Access System (UK Rail)	\$4,995
#A-92000D	Shoulder Access System (Denyer Rail)	\$4,995
#A-92001	Shoulder Access System - Sterile Disposables (Case of 10)	\$500

