


Use the Power of Your OR Table

The Table-Powered Beach Chair easily attaches to the leg section of your surgical table and positions the patient using the power of the table. Effortlessly positions patients throughout the 0° - 90° articulation range. Patient Weight Capacity: 400 lbs (181 kg).


- **Power Lift:** Utilizes the power of the surgical table to raise and lower the chair
- **Patient Weight Capacity: 400 lbs (181 kg):** Accommodates most patients
- **Sliding Headrest:** Reduces stress on the neck when raising and lowering the patient
- **Width Extension:** Adjustable width - fits all surgical tables 21³/₈" to 23" (54.3 cm - 58.4 cm) wide
- **Head Rest Design:** Allows greater movement to accommodate shorter patients
- **Unique Table Pad Integration:** Prevents gapping between the pad and chair
- **Sliding Headrest:** Reduces stress on neck when raising and lowering the patient
- **Expanded Surgical Access:** Large cut-out for excellent surgical access
- **Articulation Range:** 0° - 90°

#A-90000 Allen Table-Powered Beach Chair

What's Required	What's Included	 Allen Table-Powered Beach Chair #A-90000	 Reusable Pad Set #A-90020 1 Included	 Reusable Security Strap #A-90003 1 Included		
	What's Required	 Beach Chair Clamp #A-90013 (US), #A-90014 (EU/UK), #A-90013D (Denyer), #A-90013JP (Jpn), 2 Required Specify Rail Type				
	Counter Traction	 Counter-Traction Strap #A-90015 1 Required	OR	 Counter-Traction Brace #A-90001 1 Required	AND	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn), #A-40025 (Den) 1 Required, Specify Rail Type
	Headrest / Disposables	 Universal Head Positioner #A-90022	AND	 Universal Head Restraint Disposables #A-90023 1 Case of 10		
	Choose an Arm Support	 LPS Arm Support #A-90002 1 Required	OR	 Multi-Task Armboard #O-ABMT 1 Required	AND	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn), #A-40025 (Den) 1 Required, Specify Rail Type
	What's Optional	 Beach Chair Cart #A-90100 (US), #A-90100E (EU), #A-90100UK (UK), #A-90100D (Den), #A-90100J (Jpn) Specify Rail Type		 Allen Arm Positioner #A-92000		



Unobstructed Lateral Access



Chair Lies Flat