

ZIMMER®

Polyethylene Liner Reference Chart

Cup OD Size		AA	A	B	C	D	E	F	G	H	J	K	L	M	N	P	Q
Mates With Metal Shell OD	HGP II	42	44	46	48	50	52	54	56	58	60	62	64	66	68	70	
	HGP I		40	42	44	46	48/50	52	54	56	58	60	62	64	66	68	70
	BIAS™		40	42	44	46	48/50	52	54	56	58	60	62	64	66	68	70
	TI-BAC™ II*				40	42	44	46	48	50	52	54	56	58			
	Liners available in 22, 26, 28, and 32 mm IDs.																

*Dimension shown is without pegs